



Life Sciences

# Welcome to Pall Corporation


**Swiss Cleanroom Concept Community Event**  
**24.10.2016**

**Continuously Improving Bioprocesses**

This presentation is the work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer.  
© 2016 Pall Corporation.

# About Pall Corporation

- Global provider of filtration and separation sciences enabling technology
- Deep customer applications knowledge
- Unmatched engineering capability
- Worldwide reach, local technical support
- Delivering the most innovative process solutions



# Serving Diverse End Markets

## Pall Life Sciences



### BioPharmaceuticals

- Process
  - Biotech
  - Biologics
  - Pharmaceuticals
- Process Monitoring
  - Pharma QC
  - Protein Characterization
- Laboratory Research
  - OEM Specialty Materials
  - Core Lab
  - OEM Diagnostics
  - Laboratory QC

### Medical

- Drug Delivery
  - Intravenous
  - Injectable
  - Ophthalmic
- Healthcare Water
  - Point of Use
  - Point of Entry
  - Inline
  - Equipment
- Patient Care
  - Ventilators
  - Gas Anesthesia
  - Laparoscopic Surgery

### Food & Beverage

- Beer
- Wine & Spirits
- Alcohol-free Beverages
- Food & Dairy
- Emerging Focus Areas
  - Food Safety
  - Laboratory QC

## Pall Industrial



### Process Technologies

- Chemicals
- Oil & Gas
- Polymers
- Refineries
- Nuclear/Fossil
- Alternative Energy
- Water
- Mining
- Auto/In-plant
- Primary Metals
- Pulp & Paper
- OEMs

### Aerospace

- Commercial Aero
- Military Aero
- Marine

### Microelectronics

- Semiconductors
- Electronic Components
- Data Storage
- Displays
- Graphic Arts

# Presence

## Key Commercial Offices

- **Asia:** Singapore
- **Americas:** Port Washington, USA
- **Europe:** Fribourg, Switzerland



## Key R&D/COE Technology Centers

- **Europe:** Portsmouth, UK  
Bad Kreuznach, Germany
- **Americas:** Westborough &  
Menlo Park, USA
- **Asia:** Shanghai, China  
Singapore  
Bangalore, India



## Key Manufacturing Centers

- Ilfracombe, UK
- Newquay, UK
- Pensacola, USA
- Puerto Rico, USA
- Bad Kreuznach, Ger.
- Shanghai, China





# Kleenpak® Presto Steril Konnektor

**Continuously Improving Bioprocesses**

This presentation is the work product of Pall Corporation and no portion of this presentation may be copied, published, performed, or redistributed without the express written authority of a Pall corporate officer.  
© 2016 Pall Corporation.

# Kleenpak Presto Sterile Connector


- The GENDERLESS Kleenpak Presto sterile connector allows sterile connections to be made in unclassified environments in three simple steps – the intuitive operation provides enhanced levels of sterility assurance
- Manufactured from Bisphenol-A (BPA) free polyethersulfone (PES), these devices can be used in a wide range of applications in upstream, downstream and formulation and filling
- BPA is found in plastic resins and is of concern because of possible health affects in adults and humans



# Presto Video



# Live Demo






# Manufacturing

- Automated assembly line with 100% inspection on each device
  - On line vision system will check for
    - Defects in membrane
    - Defects in weld
    - Presence of seal
- Images of 100% inspection are stored and are linked to batch and serial number of devices
- Manufacturing quality control can be traced and audited per device
- **The highest level of traceability for sterile connector in the industry!**



# Manufacturing

- Certificate of conformance
- Details the following
  - 100% of devices have been tested for defects in membrane and welding
  - Made from animal derived components (ADC) free materials
  - PES that is used has been certified BPA free
  - Material of construction compliant to USP <87>, USP <88> and USP <661>
  - Release tests conformance to:
    - USP <788> (Particulates)
    - USP <85> (Endotoxin)

 Pall Corporation

**Certificate of Test**  
For Kleenpak™ Presto Sterile Connector

We hereby certify that the Kleenpak™ Presto Sterile Connector

Part Number:	PSC1G07 PSC1G10 PSC1G06 PSC1G11 PSC1G08	Genderless Connector with 1/4" Hose Barb Genderless Connector with 3/8" Hose Barb Genderless Connector with 1/2" Hose Barb Genderless Connector with 5/8" Hose Barb Genderless Connector with 1/2" Mini Tri-Clover (MTC)
--------------	---	--

Lot Number: LLLLLL

Has been manufactured in a controlled environment and meets or exceeds the following specifications. Kleenpak™ Presto Sterile Connectors are not supplied sterile.

<b>Materials of Construction</b> The materials of construction for the Kleenpak™ Presto Sterile Connector include: <ul style="list-style-type: none"><li>• Body (fluid path): <del>Polyethersulfone</del> (PES)</li><li>• Flange* (fluid Path): <del>Polyethersulfone</del> (PES)</li><li>• Seal (fluid path): Silicone</li><li>• Membrane (fluid path): Hydrophobic <del>Polyethersulfone</del> (PES)</li><li>• Tab: <del>Polybutylene</del> Terephthalate (PBT)</li><li>• Cover: <del>Polybutylene</del> Terephthalate (PBT)</li><li>• Protective Cap: Thermoplastic Elastomer (TPE)</li></ul>	<b>Quality Criteria</b> Each Kleenpak™ Presto Sterile Connector is 100% inspected by a vision system that detects for defects in the membrane as well as the integrity of the weld between the membrane and the body.  Samples from this manufacturing lot underwent the following tests. The lot was released by Quality Control when it was verified that their respective criteria were met:  <b>Cleanliness</b> Samples of Kleenpak™ Presto Sterile Connector from this manufacturing lot meet with adequate safety margin the current limits under USP <788> Particulate Matter in Injections with effluent counts determined microscopically.  <b>Endotoxins</b> Samples of Kleenpak™ Presto Sterile Connector from this manufacturing lot meet with adequate safety margin, current requirements under USP <85> Bacterial Endotoxins Test as determined using the Limulus Amoebocyte Lysate (LAL) reagent with an aliquot from a soak solution.
---	---

*\*Mini Tri-Clover (MTC) variant only*

The Kleenpak™ Presto Sterile Connector fluid path components meet the requirements for biological reactivity, in vivo, under the current revision of the **United States Pharmacopeia** (USP) <88> for Class VI - 121°C Plastics, and in vitro, under USP <87> (Elution Test).

The components also meet the requirements under USP <661> Containers, Physicochemical Tests – Plastics.

The fluid path consists of material listed in Title 21 of the U.S. Code of Federal Regulations (CFR) parts 170-199.

This product does not contain materials of construction that are considered specified TSE or BSE risk materials according to current legislation and guidelines (reference European CPMP EMA/410/01 and Title 21 of the U.S. Code of Federal Regulations (CFR) Part 183.5). Contact Pall for further information regarding materials of construction.

In addition to the above tests, this product met manufacturing inspection standards and requirements for full traceability by lot number. This product is manufactured under a Quality System certified to ISO 9001. Consider only unopened, undamaged packages for use. Further information is available by contacting Pall.

# Validation Data

- Mechanical tests
  - Air leak test
  - Burst test
- Functional tests
  - Bacterial challenge tests
  - Flow versus differential pressure tests
  - Fluid compatibility testing
  - Extended operation
- Extractables
- Biological reactivity tests
- Shelf life



**Danke!**



Life Sciences